



Technical Compliance Sheet Laarmann CM300 ARK Scientific


Enquiry No. PN/ADM/PUR/061/2024 (Vol - II)/853 dated 16/09/2026 – Procurement of Advanced Composite Cutting Mill – C-MET, Pune

Sl. No.	Technical Requirement	Laarmann CM300 – Offered Specification	Compliance	Remarks / Deviation
1	Working principle: portable cutting mill processing the sample by strong cutting and shear force between rotor and cutting blades, then rotor and sieves; sample collected once smaller than the sieve aperture.	Portable cutting mill (mobile stand with lockable castor wheels). Sample is cut and sheared between rotating and fixed knives, crushed repeatedly between rotor and cutting blades and between rotor and sieve; once smaller than the sieve aperture it is collected in the receiving container.	COMPLIES	NO DEVIATION
2	Operation type: Cutting Type & Cross Beater Type.	Cutting Type & Cross Beater Type.	COMPLIES	NO DEVIATION
3	Sample breaking time: less than 1 second.	Less than 1 second.	COMPLIES	NO DEVIATION
4	Feed size: 60 x 80 mm or smaller.	60 × 80 mm or smaller (up to 80 mm depending on sample).	COMPLIES	NO DEVIATION
5	Peripheral speed: Cutting type 4–20 m/s or better; Cross beater type 3–25 m/s or better.	Cutting Type: 4–20 m/s; Cross Beater Type: 3–25 m/s.	COMPLIES	NO DEVIATION
6	Final fineness: 0.1 – 4 mm or better, controlled by selection of sieve size.	0.1–4 mm; final fineness controlled by selection of different bottom sieve sizes.	COMPLIES	NO DEVIATION
7	Speed: 500–4000 rpm (continuously adjustable).	500–4000 rpm, continuously adjustable.	COMPLIES	NO DEVIATION
8	Cutting blade: manganese steel or equivalent.	Manganese steel or equivalent.	COMPLIES	NO DEVIATION
9	Batch quantity: 100 g or higher.	100 g or higher (minimum 50 g; maximum 5000 ml).	COMPLIES	NO DEVIATION
10	Handling capacity: ≥ 30 kg/h.	230 kg/h.	COMPLIES	NO DEVIATION
11	Sieves: should accommodate 0.2, 0.5, 1.0, 2.0, 4.0, 6.0, 8.0, 10.0 and 20.0 mm or equivalent.	Full range of interchangeable bottom sieves: 0.2, 0.5, 1.0, 2.0, 4.0, 6.0, 8.0, 10.0 and 20.0 mm.	COMPLIES	NO DEVIATION
12	Collector capacity: 0.25 – 30 L or bigger.	Collector volume 5000 ml (5 L) with two side handles and plexi window.	COMPLIES	NO DEVIATION

13	Grinding chamber volume: stainless steel, 1.1740 steel, hardened steel.	Stainless steel, 1.1740 steel, hardened steel.	COMPLIES	NO DEVIATION
14	Grinding chamber material: stainless steel, heavy-metal-free.	Stainless steel / heavy-metal-free.	COMPLIES	NO DEVIATION
15	Protection: IP40.	IP40.	COMPLIES	NO DEVIATION
16	Noise: ≤ 80 dB.	< 80 dB (acoustic noise-reduction Bondal infeed hopper).	COMPLIES	NO DEVIATION
17	Rated power: 1.5 kW or more.	1.5 kW.	COMPLIES	NO DEVIATION
18	Power supply: 240 V, 50 Hz.	240 V, 50 Hz.	COMPLIES	NO DEVIATION
19	Instrument size (W×D×H): 600 × 725 × 1500 mm or less.	600 × 725 × 1500 mm or less (with stand).	COMPLIES	NO DEVIATION
20	Net weight: 120 kg or less.	90 kg with underframe (68 kg without underframe) – within 120 kg.	COMPLIES	NO DEVIATION
21	Sample feeding: fed from the hopper and falls into the grinding chamber.	Sample is fed from the infeed hopper and falls into the cutting / grinding chamber via guide chute.	COMPLIES	NO DEVIATION
22	Applications: cutting/shredding soft, hard, tough, fibrous samples – plastic, rubber, plastic-inorganic composites, circuit board, cable, electronic waste, non-metallic waste.	Cutting/shredding of soft, hard, tough and fibrous samples; plastic, rubber, plastic-inorganic composites, circuit board, cable, electronic waste, non-metallic waste, etc.	COMPLIES	NO DEVIATION
23	Control system: big-screen touch panel, human interface; speed adjustment from the control panel.	Big-screen touch panel / human interface; speed adjustment from the control panel.	COMPLIES	NO DEVIATION
24	Safety: electric safety switch to prevent starting when the door is open.	Solid door with electric safety switch – machine cannot start when the door is open.	COMPLIES	NO DEVIATION
25	Operation: easy to operate and clean; rotor and sieve easily installed, removed and replaced; small cutting chamber, easy to clean.	Easy operation and cleaning; front door opens sideways, rotor / infeed hopper can be taken out and the sieve is easily exchanged; small grinding chamber, cleaned easily.	COMPLIES	NO DEVIATION
26	Grinding: avoids cross-contamination; heavy-metal-free material option; rapid with low heat; exchangeable push-fit grinding insert.	Avoids cross-contamination; heavy-metal-free material option supplied; rapid grinding with low heat, suitable for temperature-sensitive materials; push-fit grinding insert supplied.	COMPLIES	NO DEVIATION
27	Sample collection: cyclone separator suitable for light or small sample volumes.	Cyclone separator supplied – suitable for light / small sample volumes, addresses heat transfer during grinding and allows easy cleaning.	COMPLIES	NO DEVIATION

28	Standard consumables & accessories: universal V-type feed hopper; cyclone separator; cutting edges (cutting & cross-beater); stainless steel bottom sieves 0.2–20 mm; collecting bottles 250 ml, 500 ml; grinding inserts 180° and 360°.	Universal V-type feed hopper; cyclone separator; cutting edges (Cutting Type & Cross Beater Type); stainless steel bottom sieves 0.2–20 mm; collecting bottles 250 ml and 500 ml; 180° and 360° grinding inserts.	COMPLIES	NO DEVIATION
29	Installation, commissioning and training: free, minimum 1–2 days at buyer's end; demonstration on polymer, composite and rubber samples; installation requirements (space, floor, utilities) to be stated.	Free installation, commissioning and training for minimum 1–2 days at C-MET, with demonstration on polymer composites and rubber. Installation requirements (space, floor condition, utilities) will be stated in the site preparation guide.	COMPLIES	NO DEVIATION
30	Documentation: operation and maintenance manual (hard & soft copy); site preparation guide; certificate that system is brand new; OEM authorisation certificate.	Operation & maintenance manual (hard and soft copy), site preparation guide with utilities requirement, brand-new / not refurbished certificate and OEM authorisation certificate will be supplied.	COMPLIES	NO DEVIATION
31	Conformity certification: mandatory CE; mandatory ISO; Make in India preferred.	CE certification, ISO certification and Make in India certificate enclosed.	COMPLIES	NO DEVIATION
32	Warranty: minimum one year fully comprehensive; faulty components repaired/replaced free; labour, shipment and packing charges not chargeable in warranty.	One year fully comprehensive warranty; faulty components repaired or replaced free of cost, including labour, shipment and packing & forwarding charges.	COMPLIES	NO DEVIATION
33	Delivery period: within 6 to 8 weeks.	Within 6 to 8 weeks.	COMPLIES	NO DEVIATION
34	AMC and support: willing to take comprehensive and non-comprehensive AMC; minimum 7 years of spares and service support after installation.	Willing to take comprehensive and non-comprehensive AMC after warranty; minimum 7 years of spares and service support guaranteed after installation and commissioning.	COMPLIES	NO DEVIATION
35	Testimonials and feedback: list of installations of similar machines with customer contact details.	List of installations with customer contact details enclosed; customers may be contacted for feedback.	COMPLIES	NO DEVIATION
36	Machine breakdown to be attended within 48 hours.	Breakdown attended within 48 hours.	COMPLIES	NO DEVIATION
37	Demonstration of the proposed machine, in person or online, as part of technical evaluation.	Demonstration of the proposed machine will be arranged in person or online.	COMPLIES	NO DEVIATION

For ARK Scientific,



Yours Sincerely,

Aneep. Garneni
CEO

Phone

+91 9010631144

Email

aneep@arkscientific.in

Website

www.arkscientific.in

Date:02/10/2026

NO DEVIATION CERTIFICATE

To,

The In-charge Administration,
Centre for Materials for Electronics Technology (C-MET),
Panchwati, Off. Dr. Homi Bhabha Road, Pashan,
Pune – 411 008

Enquiry No.: PN/ADM/PUR/061/2024 (Vol - II)/853 dated 16/09/2026

Subject: Procurement of Advanced Composite Cutting Mill (Qty: 01 No.)

Sir / Madam,

We, **M/s ARK SCIENTIFIC**, Dhriti Towers, 5-83, near Naina Convention Road, Prashanti Nagar, IDA Kukatpally, Hyderabad, Telangana – 500072, hereby certify that the equipment offered by us, the **Laarmann CM300 Power Cutting Mill**, conforms to all the technical specifications (Annexure – I, Sl. Nos. 1 to 37) and to all the Specific Terms and Conditions of the above enquiry, as detailed in our Technical Compliance Sheet.

There is NO DEVIATION from the technical specifications, terms and conditions of the enquiry.

We also understand that any false declaration will lead to rejection of our offer and action as per the terms of the enquiry.

For ARK Scientific,



Ark Scientific
Aneep Garneni
CEO

Yours Sincerely,

Aneep. Garneni
CEO

Phone
+91 9010631144

Email
aneep@arkscientific.in

Website
www.arkscientific.in