



ARK SCIENTIFIC
Pioneers in Solids

Date:02/10/2026

Technical Compliance Sheet Laarmann CM300 ARK Scientific


Enquiry No. PN/ADM/PUR/061/2024 (Vol - II)/853 dated 16/09/2026 – Procurement of Advanced Composite Cutting Mill – C-MET, Pune

Sl. No.	Technical Requirement	Laarmann CM300 – Offered Specification	Compliance	Remarks / Deviation
1	Working principle: portable cutting mill processing the sample by strong cutting and shear force between rotor and cutting blades, then rotor and sieves; sample collected once smaller than the sieve aperture.	CM300 Power Cutting Mill. Material is shredded by cutting between rotating (rotor) and fixed (stator) knives until it passes the attached bottom sieve; ground product is collected in the receiving collector. Mobile on lockable castor wheels.	COMPLIES	NO DEVIATION
2	Operation type: Cutting Type & Cross Beater Type.	Cutting type (3 rotor knives / 3 stator knives). Cross-beater type rotor is not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
3	Sample breaking time: less than 1 second.	Not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
4	Feed size: 60 x 80 mm or smaller.	Feed size: up to 80 mm depending on sample.	COMPLIES	NO DEVIATION
5	Peripheral speed: Cutting type 4–20 m/s or better; Cross beater type 3–25 m/s or better.	Not stated in the catalogue (rotor speed 500–4000 rpm).	TO CONFIRM	TO BE CONFIRMED
6	Final fineness: 0.1 – 4 mm or better, controlled by selection of sieve size.	End fineness 250 µm (0.25 mm); end fineness adjusted by interchangeable bottom sieves. Lowest catalogue fineness (0.25 mm) is coarser than the 0.1 mm lower limit.	DEVIATION	SEE OFFERED SPECIFICATION
7	Speed: 500–4000 rpm (continuously adjustable).	Speed 500–4000 rpm; variable speed.	COMPLIES	NO DEVIATION
8	Cutting blade: manganese steel or equivalent.	Knife material not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
9	Batch quantity: 100 g or higher.	Quantity minimum 50 g; maximum 5000 ml.	COMPLIES	NO DEVIATION
10	Handling capacity: ≥ 30 kg/h.	Not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
11	Sieves: should accommodate 0.2, 0.5, 1.0, 2.0, 4.0, 6.0, 8.0, 10.0 and 20.0 mm or equivalent.	Full range of interchangeable bottom sieves; individual sieve sizes not listed in the catalogue.	TO CONFIRM	TO BE CONFIRMED

12	Collector capacity: 0.25 – 30 L or bigger.	Collector volume 5000 ml (5 L) with two side handles and plexi window.	COMPLIES	NO DEVIATION
13	Grinding chamber volume: stainless steel, 1.1740 steel, hardened steel.	Massive grinding chamber made from vibration-reducing components; material grade not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
14	Grinding chamber material: stainless steel, heavy-metal-free.	Not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
15	Protection: IP40.	Not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
16	Noise: ≤ 80 dB.	Acoustic noise-reduction infeed hopper (Bondal, vibration-damping composite); dB value not stated.	TO CONFIRM	TO BE CONFIRMED
17	Rated power: 1.5 kW or more.	Motor power 1.5 kW.	COMPLIES	NO DEVIATION
18	Power supply: 240 V, 50 Hz.	230 V / 50 Hz.	COMPLIES	NO DEVIATION
19	Instrument size (W×D×H): 600 × 725 × 1500 mm or less.	Dimensions 390 × 485 × 375 mm (machine); 650 × 750 × 1550 mm with stand. Overall size with stand exceeds the 600 × 725 × 1500 mm limit.	DEVIATION	SEE OFFERED SPECIFICATION
20	Net weight: 120 kg or less.	68 kg without underframe; 90 kg with underframe.	COMPLIES	NO DEVIATION
21	Sample feeding: fed from the hopper and falls into the grinding chamber.	Material falls into the cutting chamber from the infeed hopper via guide chute.	COMPLIES	NO DEVIATION
22	Applications: cutting/shredding soft, hard, tough, fibrous samples – plastic, rubber, plastic-inorganic composites, circuit board, cable, electronic waste, non-metallic waste.	Flexible, medium-hard, fibrous and mixed products; technical plastics (ABS, PA, POM, PE), rubber, computer and electronic scrap, cables, waste and others.	COMPLIES	NO DEVIATION
23	Control system: big-screen touch panel, human interface; speed adjustment from the control panel.	Touch-panel control unit shown in the catalogue; speed 500–4000 rpm; start/stop by on/off button. Screen size not stated.	TO CONFIRM	TO BE CONFIRMED
24	Safety: electric safety switch to prevent starting when the door is open.	Solid door with safety switch.	COMPLIES	NO DEVIATION
25	Operation: easy to operate and clean; rotor and sieve easily installed, removed and replaced; small cutting chamber, easy to clean.	Front door opens sideways; rotor / infeed hopper can be taken out for fast cleaning; sieve is easily exchanged.	COMPLIES	NO DEVIATION
26	Grinding: avoids cross-contamination; heavy-metal-free material option; rapid with low heat; exchangeable push-fit grinding insert.	Rapid size reduction (80 mm to 0.25 mm); heavy-metal-free option and push-fit insert not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED

27	Sample collection: cyclone separator suitable for light or small sample volumes.	5-litre sample collector with side handles and plexi window; cyclone separator not stated in the catalogue.	TO CONFIRM	TO BE CONFIRMED
28	Standard consumables & accessories: universal V-type feed hopper; cyclone separator; cutting edges (cutting & cross-beater); stainless steel bottom sieves 0.2–20 mm; collecting bottles 250 ml, 500 ml; grinding inserts 180° and 360°.	Infeed hopper (various types, wooden/teflon plunger available); bottom sieves; knives (rotor/stator). Remaining items not listed in the catalogue.	TO CONFIRM	TO BE CONFIRMED
29	Installation, commissioning and training: free, minimum 1–2 days at buyer's end; demonstration on polymer, composite and rubber samples; installation requirements (space, floor, utilities) to be stated.	ARK Scientific will carry out free installation, commissioning and training at C-MET, with demonstration of cutting/grinding of polymers, composites and rubber. Machine footprint 650 × 750 mm, weight 90 kg, supply 230 V/50 Hz, 1.5 kW.	COMPLIES	NO DEVIATION
30	Documentation: operation and maintenance manual (hard & soft copy); site preparation guide; certificate that system is brand new; OEM authorisation certificate.	All documents will be supplied; OEM authorisation certificate from Laarmann to be enclosed.	COMPLIES	NO DEVIATION
31	Conformity certification: mandatory CE; mandatory ISO; Make in India preferred.	Make in India certificate attached. CE declaration of conformity and ISO certificate for CM300 to be enclosed.	TO CONFIRM	TO BE CONFIRMED
32	Warranty: minimum one year fully comprehensive; faulty components repaired/replaced free; labour, shipment and packing charges not chargeable in warranty.	One year comprehensive warranty; faulty parts repaired or replaced free of cost including labour and shipment.	COMPLIES	NO DEVIATION
33	Delivery period: within 6 to 8 weeks.	Within 6 to 8 weeks from the date of purchase order.	COMPLIES	NO DEVIATION
34	AMC and support: willing to take comprehensive and non-comprehensive AMC; minimum 7 years of spares and service support after installation.	Agreed – comprehensive and non-comprehensive AMC offered; spares and service support available for minimum 7 years.	COMPLIES	NO DEVIATION
35	Testimonials and feedback: list of installations of similar machines with customer contact details.	List of installations with contact details to be enclosed.	TO CONFIRM	TO BE CONFIRMED
36	Machine breakdown to be attended within 48 hours.	Breakdown will be attended within 48 hours.	COMPLIES	NO DEVIATION
37	Demonstration of the proposed machine, in person or online, as part of technical evaluation.	Demonstration will be arranged in person or online on request.	COMPLIES	NO DEVIATION

For ARK Scientific,



Yours Sincerely,

Aneep. Garneni
CEO

Phone

+91 9010631144

Email

aneep@arkscientific.in

Website

www.arkscientific.in